

Human
CMP-N-acetylneuraminate-beta-galactosamide-alpha-2,3-sialyltransferase 1 (ST3GAL1) ELISA Kit

Catalog No: #EK20337

Package Size: #EK20337-1 48T #EK20337-2 96T

Description	
Product Name	Human CMP-N-acetylneuraminate-beta-galactosamide-alpha-2,3-sialyltransferase 1 (ST3GAL1) ELISA Kit
Species Reactivity	Human

Application Details	
Detect Range:	0.156-10 ng/mL
Sensitivity:	0.061ng/mL
Sample Type:	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume:	50 μL
Assay Time:	1-4.5h
Detection wavelength:	450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.