

Human Platelet Factor 3 (PF3) ELISA Kit

Catalog No: #EK20338

Package Size: #EK20338-1 48T #EK20338-2 96T



Description

Product Name	Human Platelet Factor 3 (PF3) ELISA Kit
Species Reactivity	Hu

Application Details

Detect Range:0.312-20ng/mL
Sensitivity:0.113ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.