

General Resolvin D1 (RvD1) ELISA Kit

Catalog No: #EK20339

Package Size: #EK20339-1 48T #EK20339-2 96T



Description

Product Name	General Resolvin D1 (RvD1) ELISA Kit
Species Reactivity	General

Application Details

Detect Range:62.5-4000 pg/mL

Sensitivity:25.6pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.