

# Rabbit Interleukin 17 (IL17) ELISA Kit

Catalog No: #EK20343

Package Size: #EK20343-1 48T #EK20343-2 96T

## Description

|                    |                                        |
|--------------------|----------------------------------------|
| Product Name       | Rabbit Interleukin 17 (IL17) ELISA Kit |
| Applications       | ELISA                                  |
| Species Reactivity | Rabbit                                 |
| Storage            | 12 months                              |

## Application Details

Detect Range:31.2-2000pg/mL

Sensitivity:11.0pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 uL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.