

Rat Ras Homolog Gene Family, Member A (RHOA)
ELISA Kit

Catalog No: #EK20367

Package Size: #EK20367-1 48T #EK20367-2 96T



Description

Product Name	Rat Ras Homolog Gene Family, Member A (RHOA) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:0.156-10ng/mL
Sensitivity:0.055ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.