

Human Collagen Type I Alpha 1 (COL1a1) ELISA Kit

Catalog No: #EK20387

Package Size: #EK20387-1 48T #EK20387-2 96T

Description

Product Name	Human Collagen Type I Alpha 1 (COL1a1) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.133ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.