

# Human Pre-mRNA-processing factor 40 homolog A (PRPF40A) ELISA Kit

Catalog No: #EK20389

Package Size: #EK20389-1 48T #EK20389-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                                   |
|--------------------|-------------------------------------------------------------------|
| Product Name       | Human Pre-mRNA-processing factor 40 homolog A (PRPF40A) ELISA Kit |
| Applications       | ELISA                                                             |
| Species Reactivity | Human                                                             |

## Application Details

Detect Range:78.1-5000pg/mL

Sensitivity:35.3pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &amp;mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.