

Human Cyclin-dependent kinases regulatory subunit 1 (CKS1B) ELISA Kit

Catalog No: #EK20390

Package Size: #EK20390-1 48T #EK20390-2 96T

Description

Product Name	Human Cyclin-dependent kinases regulatory subunit 1 (CKS1B) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:62.5-4000pg/mL

Sensitivity:26.4pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.