

Mouse cholesterol(CH) ELISA Kit

Catalog No: #EK20401

Package Size: #EK20401-1 48T #EK20401-2 96T



Description

Product Name	Mouse cholesterol(CH) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:3.12-200μg/mL

Sensitivity:0.93μg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.