

Chicken NADH Dehydrogenase, Quinone 1 (NQO1) ELISA Kit

Catalog No: #EK20403

Package Size: #EK20403-1 48T #EK20403-2 96T

Description

Product Name	Chicken NADH Dehydrogenase, Quinone 1 (NQO1) ELISA Kit
Applications	ELISA
Species Reactivity	Chicken (Gallus)
Storage	12 months

Application Details

Detect Range: 0.31-20 ng/mL

Sensitivity: 0.14 ng/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 µL

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.