

Mouse Anti-thyroid-Peroxidase Antibody (TPO-Ab) ELISA Kit

Catalog No: #EK20422

Package Size: #EK20422-1 48T #EK20422-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Anti-thyroid-Peroxidase Antibody (TPO-Ab) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:0.312-10ng/mL

Sensitivity:0.271ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.