

Human Alpha-1-Antitrypsin (a1AT) ELISA Kit

Catalog No: #EK20423

Package Size: #EK20423-1 48T #EK20423-2 96T

Description

Product Name	Human Alpha-1-Antitrypsin (a1AT) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:31.2-2000ng/mL

Sensitivity:11.37ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.