

Rat Hexokinase 2 (HK2) ELISA Kit

Catalog No: #EK20424

Package Size: #EK20424-1 48T #EK20424-2 96T



Description

Product Name	Rat Hexokinase 2 (HK2) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse
Storage	12 months

Application Details

Detect Range:1.56-100ng/mL
Sensitivity: 0.70ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.