

Canine Elastin (ELN) ELISA Kit

Catalog No: #EK20433

Package Size: #EK20433-1 48T #EK20433-2 96T

Description

Product Name	Canine Elastin (ELN) ELISA Kit
Applications	ELISA
Species Reactivity	Canine

Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.137ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.