

Rabbit Colony Stimulating Factor 3, Granulocyte (GCSF) ELISA Kit

Catalog No: #EK20435

Package Size: #EK20435-1 48T #EK20435-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rabbit Colony Stimulating Factor 3, Granulocyte (GCSF) ELISA Kit
Applications	ELISA
Species Reactivity	Rabbit

Application Details

Detect Range:7.8-500pg/mL

Sensitivity:3.1pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.