

Sheep Interferon Gamma (IFNg) ELISA Kit

Catalog No: #EK20449

Package Size: #EK20449-1 48T #EK20449-2 96T

Description

Product Name	Sheep Interferon Gamma (IFNg) ELISA Kit
Applications	ELISA Kit
Storage	12 months

Application Details

Detect Range:3.12-200pg/mL

Sensitivity:1.13pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 uL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.