

# Human Monocyte Chemotactic Protein 1 (MCP1) ELISA Kit

Catalog No: #EK20451

Package Size: #EK20451-1 48T #EK20451-2 96T

## Description

|                    |                                                       |
|--------------------|-------------------------------------------------------|
| Product Name       | Human Monocyte Chemotactic Protein 1 (MCP1) ELISA Kit |
| Applications       | ELISA                                                 |
| Species Reactivity | Human                                                 |

## Application Details

Detect Range:15.6-1000pg/mL

Sensitivity:6.6pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.