

Bovine Collagen Type III (Col 3) ELISA Kit

Catalog No: #EK20471

Package Size: #EK20471-1 48T #EK20471-2 96T

Description

Product Name	Bovine Collagen Type III (Col 3) ELISA Kit
Applications	ELISA
Species Reactivity	Bovine

Application Details

Detect Range:6.25-400 ng/mL

Sensitivity:3.52ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.