

Mouse Circulating immune complex (CIC) ELISA Kit

Catalog No: #EK20472

Package Size: #EK20472-1 48T #EK20472-2 96T

Description

Product Name	Mouse Circulating immune complex (CIC) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range: 6.25-400 ng/mL

Sensitivity: 3.13 ng/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.