

Human Kallikrein 1 (KLK1) ELISA Kit

Catalog No: #EK20475

Package Size: #EK20475-1 48T #EK20475-2 96T



Description

Product Name	Human Kallikrein 1 (KLK1) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:9.37-600pg/mL
Sensitivity:3.1pg/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.