

# Human Gasdermin-D (GSDMD) ELISA Kit

Catalog No: #EK20482

Package Size: #EK20482-1 48T #EK20482-2 96T



## Description

|                    |                                     |
|--------------------|-------------------------------------|
| Product Name       | Human Gasdermin-D (GSDMD) ELISA Kit |
| Applications       | ELISA                               |
| Species Reactivity | Human                               |

## Application Details

|                                                                                                                  |
|------------------------------------------------------------------------------------------------------------------|
| Detect Range:0.312-20 ng/mL                                                                                      |
| Sensitivity:0.105ng/mL                                                                                           |
| Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Sample Volume: 50 &mu;L                                                                                          |
| Assay Time:1-4.5h                                                                                                |
| Detection wavelength:450 nm                                                                                      |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.