

# Human Lecithin Cholesterol Acyltransferase (LCAT) ELISA Kit

Catalog No: #EK20489

Package Size: #EK20489-1 48T #EK20489-2 96T

## Description

|                    |                                                             |
|--------------------|-------------------------------------------------------------|
| Product Name       | Human Lecithin Cholesterol Acyltransferase (LCAT) ELISA Kit |
| Applications       | ELISA                                                       |
| Species Reactivity | Human                                                       |

## Application Details

Detect Range:0.468-30ng/mL

Sensitivity:0.19ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50  $\mu$ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.