

# Human Exenatide (EXT) ELISA Kit

Catalog No: #EK20492

Package Size: #EK20492-1 48T #EK20492-2 96T

## Description

|                    |                                 |
|--------------------|---------------------------------|
| Product Name       | Human Exenatide (EXT) ELISA Kit |
| Applications       | ELISA                           |
| Species Reactivity | Human                           |

## Application Details

Detect Range:125-8000pg/mL

Sensitivity:43pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.