

Human Carbonic anhydrase-related protein 11 (CA11) ELISA Kit

Catalog No: #EK20519

Package Size: #EK20519-1 48T #EK20519-2 96T

Description

Product Name	Human Carbonic anhydrase-related protein 11 (CA11) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.78-50U/mL

Sensitivity:0.29U/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.