

# Monkey Interleukin-5 (IL5) ELISA Kit

Catalog No: #EK20524

Package Size: #EK20524-1 48T #EK20524-2 96T

## Description

|                    |                                      |
|--------------------|--------------------------------------|
| Product Name       | Monkey Interleukin-5 (IL5) ELISA Kit |
| Applications       | ELISA                                |
| Species Reactivity | Monkey                               |

## Application Details

Detect Range:6.25-400pg/mL

Sensitivity:2.58pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

**Note:** This product is for in vitro research use only and is not intended for use in humans or animals.