

Mouse Creatine Kinase MB Isoenzyme (CKMB) ELISA Kit

Catalog No: #EK20526

Package Size: #EK20526-1 48T #EK20526-2 96T

Description

Product Name	Mouse Creatine Kinase MB Isoenzyme (CKMB) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:1.25-80 mU/mL

Sensitivity:0.49mU/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.