

Monkey Immunoglobulin G (IgG) ELISA Kit

Catalog No: #EK20529

Package Size: #EK20529-1 48T #EK20529-2 96T

Description

Product Name	Monkey Immunoglobulin G (IgG) ELISA Kit
Applications	ELISA
Species Reactivity	Monkey

Application Details

Detect Range:3.12-200ug/ml

Sensitivity:1.36ug/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.