

Mouse 60S acidic ribosomal protein P2 (RPLP2) ELISA Kit

Catalog No: #EK20530

Package Size: #EK20530-1 48T #EK20530-2 96T

Description

Product Name	Mouse 60S acidic ribosomal protein P2 (RPLP2) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:15.6-1000pg/mL

Sensitivity:5.8pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.