

Mouse Long-chain-fatty-acid--CoA ligase 3 (ACSL3) ELISA Kit

Catalog No: #EK20538

Package Size: #EK20538-1 48T #EK20538-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Long-chain-fatty-acid--CoA ligase 3 (ACSL3) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:15.6-1000 pg/mL

Sensitivity:6.9pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.