

Mouse Carnitine palmitoyltransferase 2, mitochondrial (CPT2) ELISA Kit

Catalog No: #EK20545

Package Size: #EK20545-1 48T #EK20545-2 96T

Description

Product Name	Mouse Carnitine palmitoyltransferase 2, mitochondrial (CPT2) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:62.5-4000 pg/mL

Sensitivity:24.5pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.