

Human Cross Linked C-Telopeptide Of Type I Collagen (CTXI) ELISA Kit

Catalog No: #EK20546

Package Size: #EK20546-1 48T #EK20546-2 96T

Description

Product Name	Human Cross Linked C-Telopeptide Of Type I Collagen (CTXI) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.059ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.