

# Mouse Bone Morphogenetic Protein Receptor 1B (BMPR1B) ELISA Kit

Catalog No: #EK20547

Package Size: #EK20547-1 48T #EK20547-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                                 |
|--------------------|-----------------------------------------------------------------|
| Product Name       | Mouse Bone Morphogenetic Protein Receptor 1B (BMPR1B) ELISA Kit |
| Applications       | ELISA                                                           |
| Species Reactivity | Mouse                                                           |

## Application Details

Detect Range:62.5-4000pg/mL

Sensitivity:24.5pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &amp;mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.