

MouseB Lipoprotein,B aB (Lpa)B ELISAB Kit

Catalog No: #EK20555

Package Size: #EK20555-1 48T #EK20555-2 96T

Description

Product Name	MouseB Lipoprotein,B aB (Lpa)B ELISAB Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:0.6ng/ml -7.2ng/ml

Sensitivity:0.12ng/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.