

Fish Interleukin 1 Beta (IL1b) ELISA Kit

Catalog No: #EK20556

Package Size: #EK20556-1 48T #EK20556-2 96T

Description

Product Name	Fish Interleukin 1 Beta (IL1b) ELISA Kit
Applications	ELISA
Species Reactivity	Fish

Application Details

Detect Range:6pg/ml-72pg/ml

Sensitivity:1.2pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.