

Mouse Carnitine Acylcarnitine Translocase (CACT) ELISA Kit

Catalog No: #EK20559

Package Size: #EK20559-1 48T #EK20559-2 96T

Description

Product Name Mouse Carnitine Acylcarnitine Translocase (CACT) ELISA Kit

Application Details

Detect Range: 0.156-10 ng/mL

Sensitivity: 0.067 ng/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 µL

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.