

Chicken Tumor Necrosis Factor Alpha (TNFa) ELISA Kit

Catalog No: #EK20561

Package Size: #EK20561-1 48T #EK20561-2 96T

Description

Product Name	Chicken Tumor Necrosis Factor Alpha (TNFa) ELISA Kit
Applications	ELISA
Species Reactivity	Chicken

Application Details

Detect Range: 7.5pg/ml-90pg/ml

Sensitivity: 1.2pg/ml

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.