

Rat Carbohydrate Antigen 15-3 (CA15-3) ELISA Kit

Catalog No: #EK20564

Package Size: #EK20564-1 48T #EK20564-2 96T



Description

Product Name	Rat Carbohydrate Antigen 15-3 (CA15-3) ELISA Kit
Host Species	Rat
Purification	ELISA

Application Details

Detect Range:7.81-500 mIU/mL
Sensitivity:3.2mIU/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.