

Mouse Ovalbumin IgG (OVA-IgG) ELISA Kit

Catalog No: #EK20573

Package Size: #EK20573-1 48T #EK20573-2 96T



Description

Product Name	Mouse Ovalbumin IgG (OVA-IgG) ELISA Kit
Host Species	Mouse
Purification	ELISA

Application Details

Detect Range:31.2-2000 pg/mL
Sensitivity:12.8pg/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.