

Mouse Low Density Lipoprotein Receptor Related Protein 6 (LRP6) ELISA Kit

Catalog No: #EK20580

Package Size: #EK20580-1 48T #EK20580-2 96T


Description

Product Name	Mouse Low Density Lipoprotein Receptor Related Protein 6 (LRP6) ELISA Kit
Host Species	Mouse
Purification	ELISA

Application Details

Detect Range:0.156-10ng/ml
Sensitivity:0.067ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.