

Porcine Alpha-2-Macroglobulin (a2M) ELISA Kit

Catalog No: #EK20581

Package Size: #EK20581-1 48T #EK20581-2 96T


Description

Product Name	Porcine Alpha-2-Macroglobulin (a2M) ELISA Kit
Purification	ELISA
Species Reactivity	Porcine

Application Details

Detect Range:1.56-100 ng/mL
Sensitivity:0.57ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.