

Human Low Density Lipoprotein Receptor (LDLR) ELISA Kit

Catalog No: #EK20583

Package Size: #EK20583-1 48T #EK20583-2 96T



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Human Low Density Lipoprotein Receptor (LDLR) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.312-20ng/mL
Sensitivity:0.118ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.