

Rat Kidney Injury Molecule 1 (Kim1) ELISA Kit

Catalog No: #EK20585

Package Size: #EK20585-1 48T #EK20585-2 96T

Description

Product Name	Rat Kidney Injury Molecule 1 (Kim1) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range: 62.5-4000 pg/mL

Sensitivity: 22.6 pg/mL


Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.