

Rat Nuclear factor- κ B p65 (NF- κ B-p65) ELISA Kit

Catalog No: #EK20587

Package Size: #EK20587-1 48T #EK20587-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Nuclear factor- κ B p65 (NF- κ B-p65) ELISA Kit
Clone No.	ELISA
Isotype	Rat

Application Details

Detect Range:62.5-4000 pg/mL

Sensitivity:25.6pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.