

Rat Alpha-1-Acid Glycoprotein (a1AGP) ELISA Kit

Catalog No: #EK20594

Package Size: #EK20594-1 48T #EK20594-2 96T

Description

Product Name	Rat Alpha-1-Acid Glycoprotein (a1AGP) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:3.12-200ng/mL

Sensitivity:1.36ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.