

Human Kallikrein 14 (KLK14) ELISA Kit

Catalog No: #EK20595

Package Size: #EK20595-1 48T #EK20595-2 96T

Description

Product Name	Human Kallikrein 14 (KLK14) ELISA Kit
Purification	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.044ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.