

Mouse N6-adenosine-methyltransferase 70 kDa subunit (METTL3) ELISA Kit

Catalog No: #EK20598

Package Size: #EK20598-1 48T #EK20598-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse N6-adenosine-methyltransferase 70 kDa subunit (METTL3) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.064ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.