

Human Interleukin 6 (IL6) ELISA Kit

Catalog No: #EK20603

Package Size: #EK20603-1 48T #EK20603-2 96T

Description

Product Name	Human Interleukin 6 (IL6) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:7.81-500pg/mL

Sensitivity:3.11pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.