

Rat TGF Beta Inducible Early Response Gene 1 (TIEG1) ELISA Kit

Catalog No: #EK20621

Package Size: #EK20621-1 48T #EK20621-2 96T

Description

Product Name	Rat TGF Beta Inducible Early Response Gene 1 (TIEG1) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.126ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.