

Mouse Glia-derived nexin(SERPINE2) ELISA kit

Catalog No: #EK20628

Package Size: #EK20628-1 48T #EK20628-2 96T

Description

Product Name	Mouse Glia-derived nexin(SERPINE2) ELISA kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:0.312-20ng/ml

Sensitivity:0.112ng/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.