

# Human Soluble cluster of differentiation 44 variant 9 (sCD44-v9) ELISA Kit

Catalog No: #EK20638

Package Size: #EK20638-1 48T #EK20638-2 96T

## Description

|                    |                                                                            |
|--------------------|----------------------------------------------------------------------------|
| Product Name       | Human Soluble cluster of differentiation 44 variant 9 (sCD44-v9) ELISA Kit |
| Applications       | ELISA                                                                      |
| Species Reactivity | Human                                                                      |

## Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.132ng/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.